

Command:AFADJ

MS2000 and RM2000 Syntax

Shortcut	AFADJ
Format	AFADJ [X=zero pot value] [Y=video amplitude value] [Z=value]
Units	integer
Remembered	Using SS Z

Tiger Syntax

Shortcut	AFADJ
Format	[Addr#]AFADJ [X=zero pot value] [Y=video amplitude value] [Z=value]
Units	integer
Type	Card-Addressed
Remembered	Using [Addr#]SS Z

X & **Y** values range between 0 and 100. Y determines amplitude of the video signal entering the system. Setting a 0 value attenuates a video signal completely, while a setting of 100 lets the full signal go through. Attenuating a video signal also reduces the noise in the signal. If the focus value on the LCD reads 2047, then the system is saturated with too much signal, try reducing the **Y** value.

X is the zeroing potentiometer; the value of X should be set such that, when **Y** is 0, the *focus value* is also 0000 (or as close as possible).

Z sets the gain of the final Analog to Digital Converter (ADC) in the auto-focus system. Range: 0 to 3. By adjusting this setting, a *focus value* for a sparse sample can be magnified to get better focus. If an incorrect value of gain is used, the ADC saturates and the focus value reaches 2047. Upon system restart, the setting returns to its default value of 0. Perform an **SS Z command** to save the current gain setting in non-volatile memory.

AFADJ	GAIN
Z=0	1x
Z=1	2x
Z=2	4x
Z=3	8x

Response

:A or Error Reply

Example

```
AFADJ X=15 Y=95
:A<CR><LF>
```

```
AFADJ
```

```
:N-3
```

(Error indicates missing arguments)

```
AFADJ X=1000 Y=-12 Z=4  
:N-4
```

(Error indicates arguments out of range)

```
AFADJ X=20 Y=90  
:N-5
```

(Error indicates operation failed, try entering one argument at a time)

```
AFADJ X? Y?  
:A X=20 Y=95
```

[commands](#), [ms2000](#), [autofocus](#)

From:

<http://www.asiimaging.com/docs/> - **Applied Scientific Instrumentation**

Permanent link:

<http://www.asiimaging.com/docs/commands/afadj>

Last update: **2024/03/15 17:57**

