

# Two-axis Control with Command Knob

More than one axis can be controlled with the command knob on the controller. Commonly, this is because in addition to a mechanical Z-axis drive, there may also be a Piezo-Z system or a Zoom axis configured in the controller.

The Liquid Crystal Display (LCD) screen shows the current position coordinates of the axes with status information displayed to the right. The W character shows which axis is controlled by the command wheel.

```
Z:  1.1234 mm :  D
F:  0.0234 mm :  W
-----
      MS2000 v8.8a
```

The front rocker switch that engages the Z-drive clutch is also used to switch the axis controlled by the command knob. The LCD screen illustrated above shows the mechanical Z-drive clutch disengaged, and the command knob controlling the F-axis. Switching the rocker switch to the other position shows:

```
Z:  1.1234 mm :  EW
F:  0.0234 mm :
-----
      MS2000 v8.8a
```

Where now the Z-drive is engaged and under control by the control knob.

[ms2000](#), [tech note](#)

From:

<http://www.asiimaging.com/docs/> - **Applied Scientific Instrumentation**

Permanent link:

[http://www.asiimaging.com/docs/two\\_axis\\_control\\_and\\_knob](http://www.asiimaging.com/docs/two_axis_control_and_knob)

Last update: **2021/09/23 17:15**

