


## **PZU-2000 Series XYZ Automated Stage with Piezo Z-axis Top Plate for Upright Microscopes**

The PZU-2000 XYZ stage has been specifically designed to provide a high resolution, and highly repeatable, means of controlling the X, Y, and Z position of the microscope stage. The XY axes derive their precise control through the use of closed-loop DC servomotors employing high-resolution rotary encoders for positioning feedback. By using closed-loop control for the stage position, there is no chance that the stage will become lost, as can occur with open-loop micro-stepped stages after a number of



moves and direction changes. The XY stage utilizes crossed-roller slides, high-precision lead screws, and zero-backlash miniature geared DC servomotors for smooth and accurate motion. The top plate of the stage accepts standard K-size slide inserts that are available for any sample, i.e., slides, petrie dishes, multi-well plates, etc. The slide insert is moved in the Z-axis via a piezo element with a range of 150  $\mu\text{m}$  with nanometer accuracy (300  $\mu\text{m}$  & 500  $\mu\text{m}$  range is also available). By moving the sample in the Z-plane, any objective can be used, eliminating twisting wires or needed spacers as required when a piezo element is put onto a single objective. The microprocessor-controlled MS-2000 control unit provides for RS-232 and USB communication with a host computer for control of the XYZ axis.

### **Features**

- Closed-loop control of the X, Y, and Z-axes for precise positioning and highly repeatable focusing
- Wide dynamic speed range with adjustable trapezoidal move profiles
- Smooth adjustable dual-range joystick control
- Backlit LCD display shows X, Y, and Z coordinates
- "Zero" and "Home" button for simple stand-alone operations
- Compact ergonomic tabletop control unit size is 6"D x 9"W x 3"H (9 x 23 x 16½ cm)
- Proven operation with many popular software packages

***We Create Solutions***

Applied Scientific Instrumentation, Inc. ♦ 29391 W Enid Rd ♦ Eugene, OR USA 97402-9533  
(541) 461-8181 ♦ (800) 706-2284 ♦ [info@ASImaging.com](mailto:info@ASImaging.com) ♦ [www.ASImaging.com](http://www.ASImaging.com)

## XY Specifications for Standard Configuration (with 6.35 mm pitch lead screws on XY stage)

XY axis range of travel	120 mm x 110 mm
XY axis resolution (encoder step)	0.088 $\mu\text{m}$
XY axis lead screw accuracy	0.25 $\mu\text{m}/\text{mm}$
XY axis RMS repeatability	< 0.7 $\mu\text{m}$
XY axis maximum velocity	7 mm/sec



## ADEPT Piezo Controller Specifications

Specification	PZU-2150	PZU-2300	PZU-2500
Piezo Travel Range (+/- 5%)	150 $\mu\text{m}$	300 $\mu\text{m}$	500 $\mu\text{m}$
Piezo smallest move / resolution*	2.2 nm	4.5 nm	7.6 nm
Maximum Load for full range travel	2 Kg	1 Kg	1 Kg
Transient Response time (10%-90%) for moves below 30% travel range**	11 – 15 ms		
External Analog input (BNC)	0 to 10 Volts		
Maximum Input Frequency	20 Hz		
Maximum Continuous Output Current	13 mA		

\*\* Time taken to travel 10%-90% for moves below 30% travel range with 600 grams load.

\* In external input mode, use of a higher bit DAC will increase resolution. For example a 0-10 analog voltage from the DAC results in the following:

PZU-2150		
External Analog Input	Steps	Resolution
16 Bit DAC	65536	2.2 nm
17 Bit DAC	131075	1.1 nm
18 Bit DAC	262144	0.55 nm

## PZU-2000 Options

- XY axes Linear Encoders for high-accuracy positioning. Linear encoder resolution is 10 nm, with a scale accuracy of 0.3  $\mu\text{m}$  per 10mm and 3  $\mu\text{m}$  per 100mm. Positioning resolution at sample is < 50 nm.
- Auto Focus (requires NTSC or PAL composite video signal).
- ASI's proven line of Z-axis drives can also be added to the fine focus shaft of the microscope to provide Z-axis positioning with a resolution of 50 nm throughout the range of the microscope's travel. The piezo unit can then be used for fast and accurate Z-axis positioning to any point within the range of travel.
- Other lead screw pitches are available for faster XY translation, or for more precise positioning when using standard rotary encoders.

***We Create Solutions***

Applied Scientific Instrumentation, Inc. ♦ 29391 W Enid Rd ♦ Eugene, OR USA 97402-9533  
(541) 461-8181 ♦ (800) 706-2284 ♦ info@ASImaging.com ♦ www.ASImaging.com