

Atmospheric Specifications for ASI Tiger Controllers

Operating	
Temperature	10° to 35° celsius
Humidity	0-75% RH non condensing
Atmospheric Pressure	700 HPa to 1060 HPa
Storage	
Temperature	-20° to 45° celsius †
Humidity	10-90% RH non condensing
Atmospheric Pressure	500 HPa to 1060 HPa

(†) must let product return to room temperature over 24 hour period before operating.

Note1: The above specs are for the ASI Tiger Controllers. Other ASI parts may have different Specifications - please inquire for individual part specs.

Note2: These are calculated specs going on ratings from individual components utilized. ASI does not test extreme conditions.

[safety](#)

From:

<https://www.asiimaging.com/docs/> - **Applied Scientific Instrumentation**

Permanent link:

https://www.asiimaging.com/docs/atmospheric_spec

Last update: **2021/09/23 17:15**

