

Command:ARM

Shortcut	ARM
Format	ARM
Remembered	Using SS Z

Sending the ARM command without arguments starts the self-scanning of the ARRAY_MODULE sequence.

This has the same effect as the [ARRAY](#) command without arguments. It is recommended that you use the ARRAY command instead of the ARM command.

This command was previously related to the deprecated [sequencer](#) module.

[commands](#), [ms2000](#)

From:

<https://asiimaging.com/docs/> - **Applied Scientific Instrumentation**

Permanent link:

<https://asiimaging.com/docs/commands/arm>

Last update: **2023/03/28 13:31**

