

# Command:ZF

Not implemented. This command was known as ZFLOCK (ZF). Command not present on MS2000 and unused on Tiger. Instead we use SLAVE\_F\_T0\_Z module that doesn't have the ability to be controlled via a command, though if needed the ability to turn on and off or change which axis is the master then it could be added.

[deprecated](#), [zflock](#)

From:

<https://www.asiimaging.com/docs/> - **Applied Scientific Instrumentation**

Permanent link:

<https://www.asiimaging.com/docs/commands/zflock>

Last update: **2024/06/25 05:02**

