

FTP Controller



The toggle switch near the joystick provides for a rapid means of moving the stage up & down. Fine control of z position is accomplished with the Z focus knob, or by computer control.



There is an option to control the Z focus with foot switches as well. If this option is requested the 9 pin serial out connector is used to connect the foot switches, and is labeled as such.

Please refer to [ZF LOCK \(FTP firmware\) module Operation](#) page for more info on FTP controller firmware.

[manual](#), [ms2000](#), [ftp](#)

From:

<http://www.asiimaging.com/docs/> - **Applied Scientific Instrumentation**

Permanent link:

http://www.asiimaging.com/docs/ftp_controller

Last update: **2021/09/23 17:15**

