

Multi-axis function

The multi-axis firmware module lets the controller generate patterns like a circle and spiral. It is implemented but has not yet undergone extensive testing. It can also be used for XY stages if the firmware module is present.

- [Command: MULTIMV \(MM\)](#) 2016/03/16 14:36
[tiger](#), [ms2000](#), [advanced feature](#)

From:

<https://asiimaging.com/docs/> - **Applied Scientific Instrumentation**

Permanent link:

<https://asiimaging.com/docs/multiaxis>

Last update: **2025/08/27 09:57**

