





# C60-TUBE-250



ITEM NO.	PART NUMBER	QTY.
1	C60-TUBE_E LENS HOLDER	1
2	C60-TUBE_E MOUNT	1
3	C60-TUBE_E LENS RETAINER	1
4	C60-TUBE-100-SPACER ONLY USED IN THE 200mm and 250mm TUBE ASSEMBLIES	1
5	C60-TUBE_E-400	1



 <b>PROPRIETARY AND CONFIDENTIAL</b> THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF APPLIED SCIENTIFIC INSTRUMENTATION. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF APPLIED SCIENTIFIC INSTRUMENTATION IS PROHIBITED.	<b>APPLIED SCIENTIFIC INSTRUMENTATION</b> 29391 W ENID RD EUGENE, OR 97402-9533 PHONE: (541)461-8181 FAX: (541)461-4018		UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES .xx ± .01 .xxx ± .005 .xxxx ± .0002 FRACTIONS ± 1/32 ANGLES ± 1°		<b>DESCRIPTION</b>  <h2>TUBE-250 ASSEMBLY</h2>	
	<b>DRAWN BY/DATE</b> M.ANTHONY 2-10-20	<b>CHECKED BY/DATE</b> -	<b>SCALE</b> 1:2	<b>REV</b> -	<b>MATERIAL</b> -	<b>DRAWING #</b> C60-TUBE-250

5

4

3

2

1