


CAM-SUP-K7-xx

(xx FOR RISER H)

OPTIONS WHEN ORDERING:

SPECIFY HEIGHT OF RISER BAR. 13", 14", 17", 19", OR 23"
 DEFAULT IS A 17" RISER BAR



PROPRIETARY AND CONFIDENTIAL
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF APPLIED SCIENTIFIC INSTRUMENTATION. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF APPLIED SCIENTIFIC INSTRUMENTATION IS PROHIBITED.

APPLIED SCIENTIFIC INSTRUMENTATION
 29391 W ENID RD
 EUGENE, OR 97402-9533
 PHONE: (541)461-8181
 FAX: (541)461-4018

UNLESS OTHERWISE SPECIFIED
 DIMENSIONS ARE IN INCHES
 .xx ± .01
 .xxx ± .005
 .xxxx ± .002
 FRACTIONS ± 1/32
 ANGLES ± 1°

DESCRIPTION

CAMERA KIT #7

POST MOUNT KIT
 WITH METRIC HOLE PATTERN CONNECTED TO BAR.

DRAWN BY/DATE

M.ANTHONY 10/15/2021

CHECKED BY/DATE

V.MARSO 10/15/2021
 J.DANIELS 10/15/2021

SCALE

1:4

REV

-

MATERIAL

303 S.S
 6061, BLK ANODIZE

DRAWING #

CAM-SUP-K7-xx