



FLAT TOP STAGE FOR OLYMPUS IX 70/73/83



PART #	DESCRIPTION
S319-2200B	Standard 6.35mm pitch leadscrews
S319-2200A	16 TPI pitch leadscrews
S319-2230A	Scan-optimized (16 TPI and no anti-backlash on X axis)



PROPRIETARY AND CONFIDENTIAL
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF APPLIED SCIENTIFIC INSTRUMENTATION. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF APPLIED SCIENTIFIC INSTRUMENTATION IS PROHIBITED.

APPLIED SCIENTIFIC INSTRUMENTATION
 29391 W ENID RD
 EUGENE, OR 97402-9533
 PHONE: (541)461-8181
 FAX: (541)461-4018

UNLESS OTHERWISE SPECIFIED
 DIMENSIONS ARE IN INCHES
 .xx ± .01
 .xxx ± .005
 .xxxx ± .0002
 FRACTIONS ± 1/32
 ANGLES ± 1°

DESCRIPTION
FLAT TOP STAGE
FITS OLYMPUS IX 70/73/83

DRAWN BY/DATE
 M.ANTHONY 2-22-21

CHECKED BY/DATE
 V.MARSO
 J.DANIELS 2-22-21

SCALE
 1:4

REV
 -

MATERIAL
 -

DRAWING #
 S319-2200B