


MIM3-PSL-O2Sxxx-2PIII

MIM3 w/ automated 2 position objective slide, two C60-CUBE-III's.
 Provides >15mm focus travel. DV6010 dovetail mount,
 Also includes top C60-SHORT-PORT-Lx, frequently used for CRISP coupling.

PART #	DESCRIPTION
O2SRMS	with RMS (Olympus) objective adapters
O2SM25	with M25 x 0.75 Threads (Nikon) objective adapters
O2SM26	with M26 x 32 TPI Threads (Mitutoyo) objective adapters
O2SM27	with M27 x 0.75 TPI Threads (Zeiss) objective adapters



PROPRIETARY AND CONFIDENTIAL
 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF APPLIED SCIENTIFIC INSTRUMENTATION. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF APPLIED SCIENTIFIC INSTRUMENTATION IS PROHIBITED.

APPLIED SCIENTIFIC INSTRUMENTATION
 29391 W ENID RD
 EUGENE, OR 97402-9533
 PHONE: (541)461-8181
 FAX: (541)461-4018

UNLESS OTHERWISE SPECIFIED
 DIMENSIONS ARE IN INCHES
 .xx ± .01
 .xxx ± .005
 .xxxx ± .0002
 FRACTIONS ± 1/32
 ANGLES ± 1°

DESCRIPTION
 MIM-CUBE-based Microscope with upper Short-Port

DRAWN BY/DATE	CHECKED BY/DATE	SCALE	REV	MATERIAL	DRAWING #
M.ANTHONY 4-1-21	J.DANIELS V.MARSO 4-1-21	1:4	0	VARIES	MIM3-PSL-O2Sxxx-2PIII